

Product Datasheet

Anti-h GFAP R13505 SPTN-5  
140060

Product Name	Anti-h GFAP R13505 SPTN-5
Catalog Number	140060
Description	Recombinant rabbit antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	Rabbit IgG
Clonality	Monoclonal
Epitope	Between amino acids 60-383 of human GFAP.
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	5.9 – 6.7
Cross Reactivity	Does not recognize vimentin, desmin or peripherin.
Specificity	Antibody recognizes human GFAP
Matched Pairs	Following pairs are especially recommended for the below mentioned assays: FIA: R13501 (capture) – R13505 (detection), and R13505 – R13504.
Shelf Life	Unspecified