

Product Datasheet

Anti-h PIGF R13204 SPTNP-5 140083

Product Name	Anti-h PIGF R13204 SPTNP-5
Catalog Number	140083
Description	Recombinant chimeric antibody containing variable region from mouse IgG and constant region from human IgG1, cultured in vitro under conditions free from animal-derived components.
Tested Applications	ELISA, CLIA, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/-10 %) (A280 nm, 1 mg/ml, 1 cm = 1.5)
Storage	+2-8°C
Isotype	Human IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.1 % ProClin 300 as a preservative
IEF Profile	7.6-9.2
Specificity	Antibody recognizes human PIGF (Placental growth factor)
Matched Pairs	Following pair is recommended for the detection of PIGF: CLIA: R13204 (capture) - R13202 or R13203 (detection)
Shelf Life	Unspecified