

Product Datasheet

Anti-h CD253 14402 SPTN-5 101024

Product Name	Anti-h CD253 14402 SPTN-5
Catalog Number	101024
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	ELISA, LF
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/-10 %)
Storage	+2-8°C
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	5.8-6.6
Specificity	Antibody recognizes the TRAIL, APO2L antigen, a 32 kDa protein
Matched Pairs	Following pairs are especially recommended for the below mentioned assays: ELISA: 14402 (capture) - 14401 (detection) LF: 14401 (membrane) - 14402 (particles)
Shelf Life	Unspecified