

## Product Datasheet

### Anti-h CD63 TS63 SPTN-5 101022

<b>Product Name</b>	Anti-h CD63 TS63 SPTN-5
<b>Catalog Number</b>	101022
<b>Description</b>	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Storage</b>	+2-8°C
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	6.5-7.6
<b>Specificity</b>	Antibody recognizes the LIMP antigen, a 30-60 kDa (smear) protein
<b>Shelf Life</b>	Unspecified