

## Product Datasheet

### Rheumatoid Factor (RF), Recombinant 610185

<b>Product Name</b>	Rheumatoid Factor (RF), Recombinant
<b>Catalog Number</b>	610185
<b>Description</b>	Recombinant rheumatoid factor (RF) monoclonal IgM expressed in hexameric form. Suitable for use as calibrator or control material, as well as in interference assays for autoimmune diagnostics.
<b>Species/Host</b>	CHO Cells
<b>Tested Applications</b>	Control, Calibration, Interference Studies
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	≥ 2,000IU/ml
<b>Storage</b>	+2-8 °C
<b>Buffer</b>	20 mM sodium phosphate, 0.5 M NaCl, pH 7.2, 0.095 % NaN <sub>3</sub> as a preservative
<b>Shelf Life</b>	12 months