

Product Datasheet

Rheumatoid Factor (RF), Recombinant 610189

Product Name	Rheumatoid Factor (RF), Recombinant
Catalog Number	610189
Description	Recombinant rheumatoid factor (RF) monoclonal IgM expressed in pentameric form. Suitable for use as calibrator or control material, as well as in interference assays for autoimmune diagnostics.
Species/Host	CHO Cells
Tested Applications	Control, Calibration, Interference Studies
Brand	Medix Biochemica
Form/Appearance	Frozen liquid
Concentration	$\geq 4,000\text{IU/ml}$
Storage	-20 °C
Buffer	20 mM sodium phosphate, 0.5 M NaCl, pH 7.2, 0.095 % NaN ₃ as a preservative
Shelf Life	12 months