

## Product Datasheet

### Anti-h CD138 13902 SPTN-5 101017

<b>Product Name</b>	Anti-h CD138 13902 SPTN-5
<b>Catalog Number</b>	101017
<b>Description</b>	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
<b>Tested Applications</b>	LF, ELISA, CLIA
<b>Immunogen</b>	U266 cell line
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	5.0 mg/ml (+/-10 %)
<b>Storage</b>	+2-8°C
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 %
<b>Affinity constant</b>	KA = 6.2 x 10 <sup>8</sup> 1/M; KD = 2.2 x 10 <sup>-9</sup> M (= 2.2 nM)
<b>Buffer</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	6.2-7.0
<b>Specificity</b>	Antibody recognizes Syndecan-1, the heparan sulfate proteoglycan, a 65-70 kDa protein
<b>Matched Pairs</b>	Following pairs are especially recommended for the below mentioned assays: ELISA: 13901 (capture) - 13902 (detection) LF: 13902 (membrane) - 13901 (particles), and 13901 - 13902 CLIA: 13902 (magnetic particles) - 13901 (detection)