

## Product Datasheet

### LH Antigen, Recombinant 610177

<b>Product Name</b>	LH Antigen, Recombinant
<b>Catalog Number</b>	610177
<b>Description</b>	Recombinant human luteinizing hormone (LH) protein, including both alpha and beta subunits, with a C-terminal histidine tag in the beta subunit.
<b>Species/Host</b>	CHO Cells
<b>Tested Applications</b>	Control, Calibration
<b>Alternative Names</b>	Luteinizing hormone, LH, Lutropin, Interstitial Cell-Stimulating Hormone
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Frozen liquid
<b>Storage</b>	-20 °C
<b>Associated Products</b>	Anti-h LH 5301, 5302, 5303, 5304 and Anti-h Alpha subunit 5501 and 6601
<b>Buffer</b>	50 mM Tris-HCl, 150 mM NaCl, pH 7.5
<b>Molecular Weight</b>	Predicted molecular weight: 24 kDa.
<b>Shelf Life</b>	12 months