

Product Datasheet

Anti-Human CD34 Class III Monoclonal Antibody, Azide Free Clone B-F23 854.260.005

Product Name	Anti-Human CD34 Class III Monoclonal Antibody, Azide Free Clone B-F23
Catalog Number	854.260.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IHC, FC
Immunogen	KG-1a cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Gp105-120 antigen, a 105-120 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells