

## Product Datasheet

### Hemoglobin A1c (HbA1c) Antigen 325-08

**Product Name** Hemoglobin A1c (HbA1c) Antigen

**Catalog Number** 325-08

**Tested Applications** Control, Calibration

**Brand** Lee Biosolutions

**Form/Appearance** Clear, red, frozen solution in phosphate buffered saline, pH 7.4

**Concentration**  $\geq 1$  mg/mL Hb (E0.1%415nm = 7)

**Storage** -20 °C

**Purity**  $\geq 96\%$  HbA1c (HPLC)

**Source** Human red blood cells

**Molecular Weight** ~65 kDa

**Shelf Life** 12 months