

Product Datasheet

Hemoglobin A1c (HbA1c), High A1c 325-21

| | |
|----------------------------|---|
| Product Name | Hemoglobin A1c (HbA1c), High A1c |
| Catalog Number | 325-21 |
| Tested Applications | Control, Calibration |
| Brand | Lee Biosolutions |
| Form/Appearance | Clear, red, frozen solution in phosphate buffered saline, pH 7.4 |
| Concentration | Total Hb: 12.0 - 17.0 g/dL (E0.1%415nm = 7) A1c %: 13.0 - 17.0% (Siemens Dimension) |
| Storage | -20 °C |
| Source | Human red blood cells |
| Molecular Weight | ~65 kDa |
| Shelf Life | 12 months |