

Product Datasheet

Hemoglobin A1c (HbA1c), High A1c 325-21

Product Name	Hemoglobin A1c (HbA1c), High A1c
Catalog Number	325-21
Description	0.2µm filtered, frozen solution in phosphate buffered saline, pH 7.4
Tested Applications	Control, Calibration
Brand	Lee Biosolutions
Form/Appearance	Clear, red solution
Concentration	Total Hb: 12.0 - 17.0 g/dL (E0.1%415nm = 7) A1c %: 13.0 - 17.0% (Siemens Dimension)
Storage	-20 °C
Source	Human red blood cells
Molecular Weight	~65 kDa
Shelf Life	12 months