

Product Datasheet

Xylazine-BSA Antigen 170106

Product Name	Xylazine-BSA Antigen
Catalog Number	170106
Description	Synthetic Xylazine, BSA conjugated
Tested Applications	LF, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-BSA conjugate
Concentration	Lot-specific
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Associated Products	140080
Buffer	Purified water
Matched Pairs	140080
Shelf Life	5 years