

Product Datasheet

Anti-Human CD34 Class III Monoclonal Antibody, PE Conjugated Clone B-F23 854.262.010

Product Name	Anti-Human CD34 Class III Monoclonal Antibody, PE Conjugated Clone B-F23
Catalog Number	854.262.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	KG-1a cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5 % BSA and 0.1 % sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Gp105-120 antigen, a 105-120 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells