

Product Datasheet

Anti-Human CD154 Monoclonal Antibody, Azide Free Clone B-B29 854.290.005

Product Name	Anti-Human CD154 Monoclonal Antibody, Azide Free Clone B-B29
Catalog Number	854.290.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	CD40L transfected 3T6 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the CD40L, TRAP, T-BAM antigen, a 28/30/33 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells