

Product Datasheet

Anti-Human CD154 Monoclonal Antibody, PE Conjugated Clone B-B29 854.292.010

Product Name	Anti-Human CD154 Monoclonal Antibody, PE Conjugated Clone B-B29
Catalog Number	854.292.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	CD40L transfected 3T6 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the CD40L, TRAP, T-BAM antigen, a 28/30/33 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells