

Product Datasheet

Anti-Human CD45 Monoclonal Antibody, Azide Free Clone B-A11 854.300.000

Product Name	Anti-Human CD45 Monoclonal Antibody, Azide Free Clone B-A11
Catalog Number	854.300.000
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Thymus cells and Jurkat cell line
Specificity	Recognises the Leucocyte Common Antigen (LCA) with different isoforms of 180-190-205-220 kDa
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells
Immunisation	Thymus cells and Jurkat cell line