

Product Datasheet

Anti-Human CD45 Monoclonal Antibody, Azide Free Clone B-A11 854.300.005

Product Name	Anti-Human CD45 Monoclonal Antibody, Azide Free Clone B-A11
Catalog Number	854.300.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	Thymus cells and Jurkat cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Leucocyte Common Antigen (LCA) with different isoforms of 180-190-205-220 kDa
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells