

Product Datasheet

Anti-Human CD45RA Monoclonal Antibody, FITC Conjugated Clone B-C15 854.321.010

Product Name	Anti-Human CD45RA Monoclonal Antibody, FITC Conjugated Clone B-C15
Catalog Number	854.321.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Reh 6 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5 % BSA and 0.1 % sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the restricted Leucocyte Common Antigen (LCA) with different isoforms of 205/220 kDa
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells