

Product Datasheet

Anti-Human CD45RO Monoclonal Antibody, Unconjugated Clone UCHL1 854.333.020

Product Name	Anti-Human CD45RO Monoclonal Antibody, Unconjugated Clone UCHL1
Catalog Number	854.333.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP, IHC
Immunogen	Human IL-2 dependent T cells
Brand	Diaclone
Form/Appearance	Phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
Storage	Store at +2-8°C. DO NOT FREEZE.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	IgG2a
Clonality	Monoclonal
Specificity	The antibody UCHL1 recognizes CD45R0, a 180 kDa low molecular weight isoform of the leukocyte common antigen (LCA). The antigen is expressed on a subset of memory/activated T cells and on cortical thymocytes.