Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 FI-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Anti-Human CD59 Monoclonal Antibody, Azide Free Clone **MEM43** 854.360

Product Name

Anti-Human CD59 Monoclonal Antibody, Azide Free Clone MEM43

854.360

Catalog Number

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications IHC, FC, IP

Immunogen Thymocytes and T lymphocytes

Brand Diaclone

Form/AppearaPhosphate-buffered saline. Sterile-filtered through 0.22 μm.

Carrier and preservative free

Stable at +2-8°C for 12 months. For longer storage freeze Storage

aliquots.

Mouse IgG2a Kappa light chain Isotype

Clonality Monoclonal

Specificity This antibody reacts with a PI-linked glycoprotein, molecular

weight 18-20 kDa, found in all types of leucocytes including

platelets.