

Product Datasheet

Anti-Human CD59 Monoclonal Antibody, Azide Free Clone MEM43 854.360

Product Name	Anti-Human CD59 Monoclonal Antibody, Azide Free Clone MEM43
Catalog Number	854.360
Description	Mouse monoclonal antibody
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IHC, FC, IP
Immunogen	Thymocytes and T lymphocytes
Alternative Names	MIC11
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	This antibody reacts with a PI-linked glycoprotein, molecular weight 18-20 kDa, found in all types of leucocytes including platelets.