

# Product Datasheet

## Anti-Human CD59 Monoclonal Antibody, Azide Free Clone MEM43 854.360.000

Product Name	Anti-Human CD59 Monoclonal Antibody, Azide Free Clone MEM43
Catalog Number	854.360.000
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IP, IHC, FC
Immunogen	Thymocytes and T lymphocytes
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	This antibody reacts with a PI-linked glycoprotein, molecular weight 18-20 kDa, found in all types of leucocytes including platelets.