

Product Datasheet

Anti-Human CD59 Monoclonal Antibody, Azide Free Clone MEM43 854.360.005

Product Name Anti-Human CD59 Monoclonal Antibody, Azide Free Clone MEM43

Catalog Number 854.360.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications IP, IHC, FC

Immunogen Thymocytes and T lymphocytes

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2a Kappa light chain

Clonality Monoclonal

Specificity This antibody reacts with a PI-linked glycoprotein, molecular weight 18-20 kDa, found in all types of leucocytes including platelets.