

# Product Datasheet

## Anti-Human CD59 Monoclonal Antibody, Unconjugated Clone MEM43 854.363.020

Product Name	Anti-Human CD59 Monoclonal Antibody, Unconjugated Clone MEM43
Catalog Number	854.363.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IP, IHC, FC
Immunogen	Thymocytes and T lymphocytes
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable for 1 month at +2-8°C. For longer storage freeze aliquots.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	This antibody reacts with a PI-linked glycoprotein, molecular weight 18-20 kDa, found in all types of leucocytes including platelets.