

Product Datasheet

Anti-Human CD82 Monoclonal Antibody, Azide Free Clone B-L2 854.440.005

Product Name	Anti-Human CD82 Monoclonal Antibody, Azide Free Clone B-L2
Catalog Number	854.440.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IHC, FC
Immunogen	CD82 transfected P815 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the CD82 antigen, a 60 kDa protein. Cross reacts with mouse CD82 protein in cytometry.
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells