

Product Datasheet

Anti-Human CD109 Monoclonal Antibody, PE Conjugated Clone B-E47 854.572.010

Product Name	Anti-Human CD109 Monoclonal Antibody, PE Conjugated Clone B-E47
Catalog Number	854.572.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	Eahy 926 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises a 170 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells