

Product Datasheet

Anti-Human CD58 Monoclonal Antibody, Azide Free Clone B-L28 854.630.005

Product Name Anti-Human CD58 Monoclonal Antibody, Azide Free Clone B-L28

Catalog Number 854.630.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications FC

Immunogen CD58 transfected CHO cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the Lymphocyte Function-associated Antigen-3 (LFA-3), a 40-70 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells