

Product Datasheet

Anti-Human CD202b Monoclonal Antibody, Azide Free Clone B-B48 854.710.000

Product Name	Anti-Human CD202b Monoclonal Antibody, Azide Free Clone B-B48
Catalog Number	854.710.000
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC, Functional assay
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Eahy 926 cell line
Specificity	Recognises Tie2, Tek, an unknown KDa protein
Biological Activity	Inhibits Tek-L / Tek binding
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells
Immunisation	Eahy 926 cell line