

Product Datasheet

Anti-Human CD14 Monoclonal Antibody, Biotin Conjugated Clone 8G3 854.774.010

Product Name	Anti-Human CD14 Monoclonal Antibody, Biotin Conjugated Clone 8G3
Catalog Number	854.774.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	Monocytes
Brand	Diaclone
Form/Appearance	Lyophilised in phosphate-buffered saline with 1% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. For longer storage freeze aliquots.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the lipopolysaccharide receptor (LPS R), a 55 kDa protein
Hybridoma	Myeloma SP2/O x Balb/c spleen cells