

Product Datasheet

Anti-Human CD261 Monoclonal Antibody, Azide Free Clone B-N36 854.850.005

Product Name	Anti-Human CD261 Monoclonal Antibody, Azide Free Clone B-N36
Catalog Number	854.850.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, Functional assay, FC
Immunogen	Recombinant human TRAIL R1/Fc chimera
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises a 43-46 kDa protein
Biological Activity	Neutralizes TRAIL induced apoptosis on HEL cell line
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells