

Product Datasheet

Anti-Human CD262 Monoclonal Antibody, Azide Free Clone B-K29 854.860.005

Product Name Anti-Human CD262 Monoclonal Antibody, Azide Free Clone B-K29

Catalog Number 854.860.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, Functional assay, FC

Immunogen Recombinant human TRAIL R2/Fc chimera

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises a 42-46 kDa protein

Biological Activity Neutralizes TRAIL induced apoptosis of U937 cell line

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells