

Product Datasheet

Anti-Human CD262 Monoclonal Antibody, Azide Free Clone B-D37 854.870

Product Name	Anti-Human CD262 Monoclonal Antibody, Azide Free Clone B-D37
Catalog Number	854.870
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC, Functional assay, IP
Immunogen	Recombinant human TRAIL R2/Fc chimera
Alternative Names	TRAIL R2 DR5
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises a 42-46 kDa protein
Biological Activity	Agonist: induced apoptosis on Jurkat cells. Enhances TRAIL induced apoptosis of U937 cell line
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells