

Product Datasheet

Anti-Human CD262 Monoclonal Antibody, Azide Free Clone B-D37 854.870.005

Product Name	Anti-Human CD262 Monoclonal Antibody, Azide Free Clone B-D37
Catalog Number	854.870.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, IP, Functional assay, FC
Immunogen	Recombinant human TRAIL R2/Fc chimera
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises a 42-46 kDa protein
Biological Activity	Agonist: induced apoptosis on Jurkat cells. Enhances TRAIL induced apoptosis of U937 cell line
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells