

Product Datasheet

Anti-Human CD263 Monoclonal Antibody, Azide Free Clone B-D44 854.890.000

Product Name Anti-Human CD263 Monoclonal Antibody, Azide Free Clone B-D44

Catalog Number 854.890.000

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications Functional assay, FC

Immunogen Recombinant human TRAIL R3/Fc chimera

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises Dcr1 (CD263), a 75.90 kDa protein

Biological Activity Inhibits TRAIL binding to TRAIL R3

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells