

Product Datasheet

Anti-Human CD263 Monoclonal Antibody, Unconjugated Clone B-D44 854.893.020

Product Name	Anti-Human CD263 Monoclonal Antibody, Unconjugated Clone B-D44
Catalog Number	854.893.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant human TRAIL R3/Fc chimera
Specificity	Recognises DcR1 (CD263), a 75.90kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells
Immunisation	Recombinant human TRAIL R3/Fc chimera