

Product Datasheet

Anti-Human Perforin Monoclonal Antibody, Azide Free
Clone B-D48
854.950

Product Name	Anti-Human Perforin Monoclonal Antibody, Azide Free Clone B-D48
Catalog Number	854.950
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	WB, IHC, ELISA, FC
Immunogen	Natural human Perforin
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognizes human Perforin, a 70 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells