

Product Datasheet

Anti-Human IL-1b Monoclonal Antibody, Azide Free Clone B-A15 855.010.005

Product Name	Anti-Human IL-1b Monoclonal Antibody, Azide Free Clone B-A15
Catalog Number	855.010.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, ELISpot, Functional assay, FC
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Store at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant antibody from the original clone
Specificity	Recognises both natural (pro and mature) and recombinant human IL-1b
Biological Activity	Inhibits IL-1β induced proliferation on D10S murine cell line
Immunisation	Recombinant human IL-1b