

## Product Datasheet

### Anti-Human IL-2 Monoclonal Antibody, Azide Free Clone B-G5 855.020

<b>Product Name</b>	Anti-Human IL-2 Monoclonal Antibody, Azide Free Clone B-G5
<b>Catalog Number</b>	855.020
<b>Description</b>	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Azide free
<b>Tested Applications</b>	WB, ELISA, FC, Functional assay, ELISpot
<b>Immunogen</b>	Natural human IL-2
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/ml (+/- 10 %)
<b>Storage</b>	Store at +2-8 °C for 12 months. DO NOT FREEZE.
<b>Isotype</b>	Mouse IgG1 Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
<b>Specificity</b>	Recognises both natural and recombinant human IL-2
<b>Biological Activity</b>	Inhibits IL-2 induced proliferation on Kit 225 cell line
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c spleen cells