

Product Datasheet

Anti-Human IL-2 Monoclonal Antibody, Azide Free Clone B-G5 855.020.005

Product Name	Anti-Human IL-2 Monoclonal Antibody, Azide Free Clone B-G5
Catalog Number	855.020.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, WB, ELISpot, Functional assay, FC
Immunogen	Natural human IL-2
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Store at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-2
Biological Activity	Inhibits IL-2 induced proliferation on Kit 225 cell line
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells