

Product Datasheet

Anti-Human IL-6 Monoclonal Antibody, PE Conjugated Clone B-E8 855.052.019

Product Name	Anti-Human IL-6 Monoclonal Antibody, PE Conjugated Clone B-E8
Catalog Number	855.052.019
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	Recombinant human IL-6
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5% BSA and 0,1% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. DO NOT FREEZE.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-6
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells