

Product Datasheet

Anti-Human IL-10 Monoclonal Antibody, Azide Free Clone B-S10 855.100.005

Product Name	Anti-Human IL-10 Monoclonal Antibody, Azide Free Clone B-S10
Catalog Number	855.100.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, Functional assay, IHC
Immunogen	Recombinant human IL-10
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-10
Biological Activity	Blocks IL-10 induced inhibition of IL-6 synthesis on LPS activated monocytes
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells