

Product Datasheet

Anti-Human IL-12 p35 Monoclonal Antibody, Azide Free Clone B-T21 855.120.005

Product Name	Anti-Human IL-12 p35 Monoclonal Antibody, Azide Free Clone B-T21
Catalog Number	855.120.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, Functional assay
Immunogen	Recombinant human IL-12
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-12 p35 + p70
Biological Activity	Inhibits IL-12 induced IFNγ synthesis on PHA activated lymphocytes
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells