

Product Datasheet

Anti-Human IL-12 p40 Monoclonal Antibody, Azide Free Clone B-P24 855.130.000

Product Name	Anti-Human IL-12 p40 Monoclonal Antibody, Azide Free Clone B-P24
Catalog Number	855.130.000
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, ELISpot, FC
Immunogen	Recombinant human IL-12
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant p40 subunit alone and in its IL-12 form (associated with p35) and in its IL-23 form (associated with p19)
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells