

Product Datasheet

Anti-Human TNF-alpha Monoclonal Antibody, Azide Free Clone B-C7 855.150.005

Product Name	Anti-Human TNF-alpha Monoclonal Antibody, Azide Free Clone B-C7
Catalog Number	855.150.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, ELISpot, Functional assay, IHC
Immunogen	Recombinant human TNF-a
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Store at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human TNF-a
Biological Activity	Inhibits TNFα induced cytotoxicity on U937 cells
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells