

Product Datasheet

Anti-Human MICA Monoclonal Antibody, Azide Free Clone B-Z26 855.370

Product Name	Anti-Human MICA Monoclonal Antibody, Azide Free Clone B-Z26
Catalog Number	855.370
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA, FC
Immunogen	Recombinant human MICA
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Store at +2-8 °C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises both natural and recombinant human MICA
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells