

Product Datasheet

Anti-Human NKp46 Monoclonal Antibody, Azide Free
Clone B-L46
855.490

Product Name	Anti-Human NKp46 Monoclonal Antibody, Azide Free Clone B-L46
Catalog Number	855.490
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c node cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Functional assay
Immunogen	Recombinant human NKp46/Fc (CD335)
Alternative Names	NCR1 CD335 Ly94
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises natural and recombinant NKp46
Biological Activity	Measured by the ability of the coated antibody to induce degranulation of activated NK cells and to enhance IFNγ secretion by non-activated NK cells
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells