

# Product Datasheet

## Anti-Human HBP Monoclonal Antibody, Azide Free Clone B-P37 855.510.005

|                     |                                                                                            |
|---------------------|--------------------------------------------------------------------------------------------|
| Product Name        | Anti-Human HBP Monoclonal Antibody, Azide Free Clone B-P37                                 |
| Catalog Number      | 855.510.005                                                                                |
| Description         | Mouse monoclonal                                                                           |
| Species/Host        | Mouse                                                                                      |
| Conjugation         | Azide free                                                                                 |
| Tested Applications | ELISA                                                                                      |
| Immunogen           | Recombinant human Azurocidin/Cap37/HBP                                                     |
| Brand               | Diaclone                                                                                   |
| Form/Appearance     | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free |
| Storage             | Stable at +2-8°C for 12 months. For longer storage freeze aliquots.                        |
| Isotype             | Mouse IgG2b Kappa light chain                                                              |
| Clonality           | Monoclonal                                                                                 |
| Specificity         | Recognises both natural and recombinant human HBP                                          |
| Hybridoma           | Myeloma X63/AG.8653 x Balb/c node cells                                                    |