

Product Datasheet

Anti-Human MCP-1 Monoclonal Antibody, Azide Free Clone B-Z21 855.580

Product Name	Anti-Human MCP-1 Monoclonal Antibody, Azide Free Clone B-Z21
Catalog Number	855.580
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c node cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human MCP-1
Alternative Names	SCYA-2 CCL2 JE
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises both natural and recombinant MCP-1 (Monocyte Chemotactic Protein1), also called CCL2 or JE
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells