

# Product Datasheet

## Anti-Human Galectin-3 Monoclonal Antibody, Azide Free Clone B-S43 855.610.005

|                     |                                                                                               |
|---------------------|-----------------------------------------------------------------------------------------------|
| Product Name        | Anti-Human Galectin-3 Monoclonal Antibody, Azide Free Clone<br>B-S43                          |
| Catalog Number      | 855.610.005                                                                                   |
| Description         | Mouse monoclonal                                                                              |
| Species/Host        | Mouse                                                                                         |
| Conjugation         | Azide free                                                                                    |
| Tested Applications | ELISA, FC                                                                                     |
| Immunogen           | Recombinant human Galectin-3                                                                  |
| Brand               | Diaclone                                                                                      |
| Form/Appearance     | Phosphate-buffered saline. Sterile-filtered through 0.22 µm.<br>Carrier and preservative free |
| Storage             | Stable at +2-8°C for 12 months. For longer storage freeze aliquots.                           |
| Isotype             | Mouse IgG2a Kappa light chain                                                                 |
| Clonality           | Monoclonal                                                                                    |
| Specificity         | Recognises both natural and recombinant protein                                               |
| Hybridoma           | Myeloma X63/AG8. x Balb/c node cells                                                          |